

## LIFEPAK® CR2 Defibrillator

### LIFEPAK CR2 Defibrillator with LIFELINKcentral™ AED Program Manager vs. ZOLL® AED 3™ Defibrillator

Physio-Control is the pioneer in defibrillation technology with more than 60 years of experience. Our complete and connected AED response system is expertly designed for user confidence. The intuitive CR2 integrates advanced technology that enables heart rhythm analysis while CPR is performed. And LIFELINKcentral AED Program Manager monitors each AED every single day, sending alerts via cellular network or Wi-Fi® about anything that may affect a unit's readiness—all automatically. This means our customers can readily overcome their most difficult challenges with a single, connected system. Here is a side-by-side comparison with the ZOLL AED 3™ AED.

ZOLL will say the Wi-Fi capabilities of the AED 3 offer seamless communication with your wireless network<sup>1</sup>, but consider having the freedom to place your AED almost anywhere in the world and be able to ensure readiness, without having to rely on a Wi-Fi network.



LIFEPAK CR2 Defibrillator with LIFELINKcentral AED Program Manager



ZOLL AED 3 (public access) AED

#### Connected. Ready.

Connection to an AED program management system via Wi-Fi or cellular network	Yes	Wi-Fi only <sup>2</sup>
Automatic tracking of AED location via the cellular network	Yes	No
Wireless monitoring of battery charge and self-tests through an AED program management system to ensure readiness	Yes	No
Email notification via Wi-Fi of device in use, including location and time	Yes	No
Wi-Fi transmission of event data to arriving emergency responders	Yes	No
Wireless distribution of AED event reports to professional caregivers	Yes (optional)	Yes <sup>3</sup>
Wireless software upgrades	Yes	No

#### Continuous high-quality CPR

Analysis during chest compressions shortens or eliminates pauses, to maximise hands-on time	Yes (with cprINSIGHT™ technology)	No
Metronome and CPR coaching to audibly guide users in performing effective compression rate and depth	Yes <sup>4</sup>	Yes
Integrated CPR feedback using only the pads, without the use of an accelerometer	Yes	No

#### Designed for user confidence

Bold and simple graphics help users focus on one task at a time	Yes	Yes
Time to shock from the end of CPR interval	Less than 7 seconds for a 360J shock (with cprINSIGHT™ technology)	Less than 10 seconds charge time for a 200J shock <sup>3</sup>
Side-by-side electrodes without a liner to minimise confusion	Yes	No
Detection of background noise and automatic volume adjustment to ensure instructions are audible	Yes	No
Optional dual language capability for use in multilingual communities	Yes	Yes
High IP rating for protection against water and dust	Yes, IP55	Yes, IP55 <sup>2</sup>
Child mode button for quick access to reduced energy	Yes	Yes
Escalating biphasic full energy for difficult-to-defibrillate patients	Yes, 150–360J	No, only escalates 120 to 200J <sup>1</sup>

## REFERENCES

- 1 ZOLL Medical Corporation. (2016). ZOLL AED 3 Spec Sheet (MCN PP 1409 0070-05).
- 2 ZOLL Medical Corporation. (2016). *ZOLL AED 3 Automated External Defibrillator Brochure* (MCN PP 1410 0068-05).
- 3 ZOLL Medical Corporation. (March 2016). *ZOLL AED 3 Administrator's Guide* (9650-000752-12 Rev. D).
- 4 Physio-Control Internal Semi-Automatic AED Comparison Usability Study, August 2016.

All claims valid as of March 2018.

If you purchased your LIFEPAK CR2 Defibrillator from an authorized Physio-Control distributor or reseller, this distributor or reseller will have access to your LIFELINKcentral AED Program Manager account and may receive notifications prompted by the LIFEPAK CR2 Defibrillator. Please note that this setting to notify your distributor or reseller can be disabled at ANY time: if you wish to disable this setting, please send a request to Physio-Control Customer Support to self-manage your site without notifications to your distributor or reseller.

To the extent that Physio-Control collects personally identifiable information through the LIFELINKcentral AED Program Manager, such information is collected pursuant to Physio-Control's privacy policy, a copy of which is located here: <http://www.physio-control.com/privacy.aspx>.

Physio-Control is now part of Stryker.

For further information, please contact Physio-Control at 800.895.5896 or visit our website at [www.physio-control.ca](http://www.physio-control.ca)

---

### Physio-Control Headquarters

11811 Willows Road NE  
Redmond, WA 98052  
[www.physio-control.com](http://www.physio-control.com)

### Customer Support

P. O. Box 97006  
Redmond, WA 98073  
Toll free 800 442 1142  
Fax 800 426 8049

### Physio-Control Canada

Physio-Control Canada Sales, Ltd.  
45 Innovation Drive  
Hamilton, ON  
L9H 7L8  
Canada  
Toll free 800 895 5896  
Fax 866 430 6115